

ThinkGear Socket Protocol

June 28, 2009



The NeuroSky product families consist of hardware and software components for simple integration of this bio-sensor technology into consumer and industrial end-applications. All products are designed and manufactured to meet exacting consumer specifications for quality, pricing, and feature sets. NeuroSky sets itself apart by providing building-block component solutions that offer friendly synergies with related and complementary technological solutions.

Reproduction in any manner whatsoever without the written permission of NeuroSky Inc. is strictly forbidden. Trademarks used in this text: eSense™, ThinkGear™, MDT™, NeuroBoy™ and NeuroSky™ are trademarks of NeuroSky Inc.

NO WARRANTIES: THE DOCUMENTATION PROVIDED IS "AS IS" WITHOUT ANY EXPRESS OR IMPLIED WARRANTY OF ANY KIND INCLUDING WARRANTIES OF MERCHANTABILITY, NONINFRINGEMENT OF INTELLECTUAL PROPERTY, INCLUDING PATENTS, COPYRIGHTS OR OTHERWISE, OR FITNESS FOR ANY PARTICULAR PURPOSE. IN NO EVENT SHALL NEUROSKY OR ITS SUPPLIERS BE LIABLE FOR ANY DAMAGES WHATSOEVER (INCLUDING, WITHOUT LIMITATION, DAMAGES FOR LOSS OF PROFITS, BUSINESS INTERRUPTION, COST OF REPLACEMENT GOODS OR LOSS OF OR DAMAGE TO INFORMATION) ARISING OUT OF THE USE OF OR INABILITY TO USE THE DOCUMENTATION PROVIDED, EVEN IF NEUROSKY HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. , SOME OF THE ABOVE LIMITATIONS MAY NOT APPLY TO YOU BECAUSE SOME JURISDICTIONS PROHIBIT THE EXCLUSION OR LIMITATION OF LIABILITY FOR CONSEQUENTIAL OR INCIDENTAL DAMAGES.

Contents

Introduction	4
Conventions	4
Packet Syntax	5
Structure	5
Header	5
Payload	5
Example Packet	6
Parsing	8
Pseudocode	8
Example Code	8
ActionScript 3.0	8

Introduction

This document defines the protocol for communicating with devices or software that implement the ThinkGear Socket Protocol (TGSP), such as the [ThinkGear Connector](#) (TGC). The TGSP was designed to allow languages and/or frameworks without a standard serial port API (e.g. Flash and most scripting languages) to integrate brainwave-sensing functionality.

TGSP streams data read from ThinkGear technology, such as the MindSet headset, to an open socket via TCP. Any language and/or framework that has a socket library will thus be able to communicate with the headset through this open socket.

Important: This document is a **draft specification** and is *extremely likely* to change prior to the final release of the document.

Conventions

Before delving into the packet specification, it is important to establish a few formatting conventions used in this document to facilitate explanation.

[TYPE]

Surrounding a TYPE with square brackets means that the TYPE is exactly a single byte long.

[TYPE. . .]

A TYPE trailed by an ellipsis indicates that the TYPE consists of one or more bytes.

([TYPE])

TYPE byte(s) surrounded by parentheses indicates that this byte may or may not be present, depending on the context.

Packet Syntax

TGSP uses a packet structure similar to that of the ThinkGear Communications Protocol (TGCP) — in fact, the `PAYLOAD` portion of the structure is largely identical. The primary difference is that TGSP does not contain a checksum of the packet, which is used in the TGCP for error-checking.

Structure

```
[ HEADER... ] [ PAYLOAD... ]
```

Every packet sent by a TGSP server will contain a fixed-length header and a variable-length payload.

Header

```
[ SYNC] [ SYNC]
```

The header for each packet consists of two `SYNC` bytes, having a value of `0xAA`. It is useful for determining the start of a packet, in the case where the implementing application loses its place in the buffer and needs to realign itself. In general, one should scan the buffer for these two `SYNC` bytes before starting to parse the payload.

Payload

```
[ CODE] ([ LENGTH]) [ DATA... ]
```

Every packet payload contains any number of the above tuple, where the exact number of tuples per packet is dictated by how the ThinkGear has been configured. The `CODE` and the optional `LENGTH` are always a single unsigned byte, whereas `DATA` can be a variable number of bytes long. The `CODE` portion of the packet dictates whether there is a `LENGTH` portion, and also how the `DATA` portion of the packet should be parsed. The following table lists the available `CODE`s, their corresponding `LENGTH` and `DATA`, and a short description of their function.

CODE	DATA content	LENGTH	Description
0x02	One (1) unsigned byte	Not used	Poor signal level
0x04	One (1) unsigned byte	Not used	eSense Attention value
0x05	One (1) unsigned byte	Not used	eSense Meditation value
0x81	Eight (8) big-endian floats	0x20	EEG powers , corresponding to delta, theta, low beta, high beta, low alpha, high alpha, low gamma, and high gamma, respectively

Note that there is **no guarantee** of payload order; EEG powers may come before the eSense values, for example.

Example Packet

Let us examine a standard packet sent by a TGSP server to an application. The raw output of the byte stream looks like the following (interpreted as hexadecimal numbers):

```
AA AA 02 1A 04 27 05 62 81 20 38 F1 50 C1 35 BD C0 55 39 0D A7
A7 38 8C 51 78 37 78 35 C6 35 3A CC CF 35 0D 61 CD 37 6C 1B 71
```

Here is the interpretation of each byte:

Index	Byte	Description
0	0xAA	SYNC
1	0xAA	SYNC
2	0x02	Poor signal CODE
3	0x1A	Poor signal level of 26
4	0x04	eSense attention CODE
5	0x27	eSense attention value of 39
6	0x05	eSense meditation CODE
7	0x62	eSense meditation value of 98
8	0x81	EEG powers CODE
9	0x20	EEG powers LENGTH
10-13	0x38F150C1	EEG delta value of 1.15×10^{-4}
14-17	0x35BDC055	EEG theta value of 1.41×10^{-6}
18-21	0x390DA7A7	EEG low alpha value of 1.35×10^{-4}
22-25	0x388C5178	EEG high alpha value of 6.69×10^{-5}
26-29	0x377835C6	EEG low beta value of 1.47×10^{-5}
30-33	0x353ACCCF	EEG high beta value of 6.95×10^{-7}
34-37	0x350D61CD	EEG low gamma value of 5.26×10^{-7}
38-41	0x376C1B71	EEG high gamma value of 1.40×10^{-5}

Parsing

By default, a TGSP server will transmit a packet of data every second (1Hz). Any parser implementations can thus keep the main parsing loop in a busy-wait state for most of the time.

Pseudocode

The pseudocode for the parser is as follows:

```
while the socket connection is open and there is data in the buffer
  scan the buffer until the sync bytes are seen

  while there is still data in the socket buffer
    read a code byte

    if the code byte refers to "attention"
      read the attention value byte
    else if the code byte refers to "meditation"
      read the meditation value byte
    else if the code byte refers to ...
      ...
    else if the code byte refers to "EEG powers"
      read eight EEG power floating point numbers
    else if the code byte is unrecognized
      continue the loop

  sleep for 1 second
```

The above pseudocode applies to most generic procedural languages. Note that in frameworks that support event handlers (such as `ActionScript`), the structure of the code can look slightly different.

Example Code

The following examples demonstrate functional implementations of TGSP parsers in a variety of languages.

ActionScript 3.0

ActionScript 3.0 (AS3) is commonly used in application frameworks from Adobe, such as Adobe Flash, Adobe Flex, and Adobe AIR. The following portion of code uses AS3's `Socket` class, which was a new introduction to ActionScript. As such, Adobe application frameworks that use older versions of ActionScript will not be able to connect to TGSS.

```
/**
 * Class constructor
```



```

*/
public function YourClass(){
    // set up a socket connection to TGSS (defaults to localhost:13854)
    thinkGearSocket = new Socket("127.0.0.1", 13854);

    // add socket event listeners
    thinkGearSocket.addEventListener(ProgressEvent.SOCKET_DATA, dataHandler);
}

/**
 * Event handler for when data is received on the socket.
 */
private function dataHandler(e : ProgressEvent){
    // buffer to store the headset data
    var buffer : ByteArray = new ByteArray();

    // read the bytes from the socket into the buffer
    thinkGearSocket.readBytes(buffer, 0, thinkGearSocket.bytesAvailable);

    buffer.position = 1; // set buffer position to anticipate SYNC bytes
    buffer.endian = Endian.BIG_ENDIAN; // set float interpretation to big-endian

    // let's look for the SYNC bytes.
    // the loop consumes the buffer until two consecutive 0xAA bytes are seen
    while(buffer.position < buffer.length &&
        buffer.readUnsignedByte() != 0xAA && buffer[buffer.position - 1] != 0xAA){
        continue;
    }

    // now process the buffer.
    while(buffer.position < buffer.length){
        // read the CODE byte
        var code : uint = buffer.readUnsignedByte();

        switch(code){
            // signal quality - read a single byte
            case 0x02:
                poorSignal = buffer.readUnsignedByte();
                break;
            // attention - read a single byte
            case 0x04:
                attention = buffer.readUnsignedByte();
                break;
            // meditation - read a single byte
            case 0x05:
                meditation = buffer.readUnsignedByte();
                break;
            // EEG powers - read the length, and eight floats (8 * 4 bytes)
            case 0x81:
                var eegPowerLength : Number = buffer.readUnsignedByte();
                delta = buffer.readFloat();
                theta = buffer.readFloat();
                beta1 = buffer.readFloat();
                beta2 = buffer.readFloat();
                alpha1 = buffer.readFloat();
                alpha2 = buffer.readFloat();
                gamma1 = buffer.readFloat();
                gamma2 = buffer.readFloat();
        }
    }
}

```

```
        break;
    }
}

// clear the buffer of any contents, in case there are several payload
// tuples that aren't handled by the parser
thinkGearSocket.flush();
}
```